
Distributed by / Vertrieb / Distribuido por / Distribué par / Distribuito da / Gedistribueerd door /  
**Дистрибутор** / Distribuije / Distribueret af / Jakelija / Διαμέτεται από / Forgalmazó / 販売元 / Dystrybutor / Distribuido por /  
 Distributor / **Дистрибутор** / Distribútor / Distributer / Distributer / Distribuerad av / Turustaja / Izplatītājs /

- |  |   |  |
|--|---|--|
| <b>GB</b> B. Braun Medical Ltd.<br>Brookdale Road<br>Sheffield, S35 2PW<br>United Kingdom                          | <b>CZ</b> B. Braun Medical s.r.o.<br>148 00 Praha 4 - Česká republika   | <b>RO</b> B. Braun Medical S.R.L.<br>Sânandrei, Str. Bernd Braun, nr. 1,<br>Jud. Timiș, 307375 - România<br>Tel. +40 256 284 905<br>contact.ro@bbraun.com<br>www.bbraun.ro |
| <b>DE</b> B. Braun Melsungen AG34212<br>Melsungen - Germany  | <b>DK</b> B. Braun Medical A/S<br>2000 Frederiksberg - Denmark  | <b>RS</b> Beohem-3 doo<br>11000 Belgrade - Serbia  |
| <b>ES</b> B. Braun Medical, S.A.U.<br>Carretera de Terrassa, 121<br>08191 Rubí (Barcelona) - España                | <b>EE</b> B.Braun Medical OÜ<br>Piltvee 8, Tallinn 12618- Estonia   | <b>RU</b> Б.Браун Медикал, 191040, г. Санкт-<br>Петербург, ул. Пушкинская, д. 10. Тел./<br>факс: (812) 320-40-04.  |
| <b>FR</b> B. BRAUN MEDICAL SAS<br>92210 Saint-Cloud - France   | <b>FI</b> B. Braun Medical Oy<br>00350 Helsinki - Finland   | <b>SK</b> B. Braun Medical s.r.o.<br>83103 Bratislava - Slovakia   |
| <b>IT</b> B. Braun Milano S.p.A.<br>20161 Milano - Italy   | <b>HU</b> B. Braun Trading Kft<br>1023 Budapest - Felhézvi utca 5<br>Magyarország   | <b>SI</b> Medis, d.o.o.<br>1001 Ljubljana - Slovenia   |
| <b>NL</b> B. Braun Medical B.V.<br>Euterpehof 10, NL-5340 AR Oss.<br>Netherlands                                   | <b>JP</b> 製造販売元 ビー・ブラウンエスクラップ株式会社<br>〒113-0033 東京都文京区本郷2-38-16  | <b>SV</b> B. Braun Medical AB<br>Box 110 - 182 12 Danderyd -<br>Sweden   |
| <b>AU</b> B. Braun Australia Pty. Ltd.<br>Level 5, 7-9 Irvine Place<br>Bella Vista NSW 2153<br>Phone: 1800 251 705 | <b>LV</b> B. Braun Medical SIA<br>Ūdelu iela 16 - Rīga, Latvija - 1064  |  |
| <b>BG</b> "Б. Браун Медикал" ЕООД<br>гр. София 1592 - България<br>тел. +359 2 807 6740                             | <b>PL</b> Aesculap-Chifa Sp. z o. o.<br>64-300 Nowy Tomyśl - Poland   |  |
|  | <b>PT</b> B. Braun Medical, Unipessoal Lda.<br>Estrada Consiglieri Pedroso 80,<br>Queluz de Baixo, 2730-053 Barcarena<br>Portugal |  |


# Flexima®/Softima® Active

5120808 / 2023-10 / M1

## Flexima®/Softima® Active

Effective

Effective


B. BRAUN MEDICAL SAS | 26 Rue Armengaud | 92210 Saint-Cloud | FRANCE

Document No.: NA - Version: 2.0 - Document ID: GMS-OMSC-000534 Print Date - Gedrukt am: 2023-12-22 19:15 (CET)

# Flexima®/Softima® Active

## EN Instructions for Use

### What Flexima® / Softima® Active is and what it is used for:

Flexima® / Softima® Active is a 1-Piece collection system for enterostomates or urostomates comprising of different types of bags:

- Closed for colostomates with formed stools,
- Drainable for colostomates with soft or liquid stools, ileostomates and jejunostomates with either "Roll'Up" outlet or "High Flow" cap (available only for Active O'Convex) allowing connection to a B. Braun distal bag,
- Drainable for urostomates with a "Uro" tap and uro connector allowing connection to a B. Braun distal bag.

The soft convex skin protector Flexima® / Softima® Active O'Convex is specifically designed to ensure a secure seal around stomas which are located in shallow to moderately deep skin folds, flush or slightly retracted.

Flexima® / Softima® Active may be used concomitantly with B. Braun stoma care accessories.

Flexima® / Softima® Active device is made of hydrocolloid skin adhesive, polymer/copolymer plastics and nonwoven film.

### Intended use:

The product is intended for collection and containment of stools or urine by patients with enterostomy or urostomy.

The product is not intended to be used on injured skin or mucus membrane.

### Indication:

The product is indicated for patients with an enterostomy and urostomy.

Flexima®/Softima® Active can only be used in adults. Not for children. Not for infants.

The intended users are healthcare professionals (HCP), caregivers and patients.

The responsible HCP must ensure that the patients and caregivers are duly trained and instructed on the use of ostomy bags, peristomal skin complications and leakage.

### Contraindications:

- Known hypersensitivity of the patient to any material of the product (see product description)
- Inability to self-stoma care in the absence of an appropriately trained caregiver or attendant for patient with poor manual dexterity and cognitive impairment.

### Side effects:

- Unknown hypersensitivity reaction to any material used leading to anaphylactic or anaphylactoid reactions.

WARNING: If severe hypersensitivity reactions occur, remove the product immediately, discontinue its use and contact a physician.

- Peristomal skin complications.

WARNING: Leakage can promote peristomal skin complications.

Note: The use of stoma care accessories to promote skin health can be considered.

### 1 Preparing the skin

Wash your hands before starting the procedure.

Gently clean the peristomal area with lukewarm water. Avoid ether, alcohol or any other topical product which could interfere with the skin protector adhesion. Carefully dab dry the skin with soft and non-shedding cloths.

Choose the appropriate size of the skin protector/ bag to protect the peristomal skin from exposure to stool and urine. Stoma size may vary directly after the surgery and at the later stage, please adapt accordingly for appropriate fitting of skin protector/bag.

### 2 Preparing the skin protector

If you are using a Cut-To-Fit product, cut the central opening of the hydrocolloid adhesive to the size of your stoma (A) just before the application, possibly with the help of the Stoma Guide. Do not cut beyond the maximum possible diameter and make sure that there is always a 2 to 3 mm margin around your stoma in order not to strangle or hurt the stomal mucosa.

Note: Pre-Cut products are not intended to be cut out.

### 3 Preparing the bag

Separate the bag walls by crumpling it slightly to introduce air into it. For drainable bag, close the closure outlet:

- Drainable "Roll'Up" bag: Fold the extremity of the sleeve four times using the guides and fold the tab back (G). To seal the outlet, make sure that the two self-fastening strips are aligned, facing each other and apply pressure on all their length.
- Drainable "High Flow" bag: Close the cap as far as maximum.
- For drainable "Uro" bag, make sure that the tap is closed by turning it clockwise. (J,K)

### 4 Applying the skin protector

Peel off the release (B). Position the bottom of the skin protector around the base of the stoma. Press the bottom part of the skin protector smoothly onto the skin, then place and press onto the upper part (the skin protector is flexible and adapts perfectly to body movements). (C)

In case of error in the positioning, it is possible to remove the bag

- EN Medical device as per 2017-745/EU
- DE Medizinprodukt gemäß 2017-745/EU
- ES Producto sanitario según el Reglamento (UE) 2017-745
- FR Dispositif médical selon 2017-745/EU
- IT Dispositivo medico secondo 2017-745/UE
- NL Medisch hulpmiddel conform 2017-745/EU
- BG Медицинско изделие съгласно 2017-745/EU
- CS Zdravotnický prostředek dle nařízení (EU) 2017/745
- DA Medicinsk udstyr i henhold til 2017-745/EU
- EL Ιατροτεχνολογικό προϊόν σύμφωνα με τον Κανονισμό 2017/745/ΕΕ
- ET Meditsiiniseade 2017-745/EU kohaselt
- FI Lääkinnällinen laite direktiivin 2017-745/EU mukaisesti
- HU Orvosi eszköz, a 2017/745/EU rendeletnek megfelelően
- JA 2017-745/EUに従った医療機器
- LV Mediciniskā ierīce saskaņā ar 2017-745/EU
- PL Wyrób medyczny zgodnie z rozporządzeniem (UE) 2017/745
- PT Dispositivo médico de acordo com a diretiva 2017-745/EU
- RO Dispozitiv medical conform 2017-745/EU
- RU Медицинское изделие согласно 2017-745/EU
- SK Zdravotnícka pomôcka podľa normy 2017-745/EU
- SL Medicinski pripomoček v skladu z 2017-745/EU
- SR Medicinsko sredstvo u skladu sa 2017-745/EU
- SV Medicinteknisk produkt enligt 2017-745/EU

# Effective

- EN Distributor
- DE Händler
- ES Distribuidor
- FR Distributeur
- IT Distributore
- NL Distributeur
- BG Дистрибутор
- CS Distributor
- DA Distributør
- EL Διανομέας
- ET Levitaja
- FI Jakelija
- HU Forgalmazó
- JA 販売代理店
- LV Izplatītājs
- PL Dystrybutor
- RO Distribuitor
- RU Дистрибутор
- SK Distribútor
- SL Distributer
- SR Distributer
- SV Distributör

# Effective

## EN Instructions for Use

and replace it properly.

In order to improve the adhesion of Flexima® / Softima® Active O'Convex products, place hands over the skin protector to warm it up.

### 5 Emptying of drainable bag

Position yourself appropriately.

#### ■ Drainable "Roll'Up" bag:

Release the two self-fastening strips and unroll the sleeve toward the toilet bowl. Be careful, emptying will start right at this point. To fully open the outlet, apply slight pressure with your fingers on guide (E). Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity (use only non-shedding material) (F).

#### ■ Drainable "High Flow" bag:

Remove the cap by pulling it straight and point it towards the toilet bowl. Once the emptying is complete, clean the cap.

#### ■ Drainable "Uro" bag

Point the tap towards the toilet bowl, slowly turn the tap counter-clockwise until the desired flow rate is obtained.

After emptying of drainable bag, close the closure outlet again as indicated above. (J,K)

### 6 Using B. Braun belt (for Flexima® / Softima® Active O'Convex only)

The soft convex skin protector helps to prevent leakages by tight-fitting the peristomal area. For extra security you can use one out of the three positions of the belt loops of the soft convex skin protector to fix the belt. Choose two loops that are diametrically opposed. (H, I)

### 7 Using filter protection stickers

Before bathing or showering, stick the filter protection sticker over the gas evacuation openings in order to maintain filter performances (D). For Split cover, this must require lifting the comfort film and to stick the filter protection filter directly on the bag films. After bath or shower, remove the filter protection stickers.

### 8 Compatibility with urostomy B. Braun distal bags

B. Braun offers a complete solution with urine collection bags for long-time intervals and during the night, such as Uribag® and Urimed® that are compatible to all B. Braun Uro bags.

The connection of the bag to a urostomy B. Braun distal bag can require the use of the supplied uro connector. To clip on the uro connector or the B. Braun Uribag®, press the two wings and insert it into the bag tap. Open the tap turning it half a turn counter-clockwise to allow urine to flow in the B. Braun distal bag. To un-clip, close the tap of the bag by turning it half clockwise, pull out the uro connector or the B. Braun Uribag® by pressing the two wings.

### 9 Compatibility with Flow Collector®

B. Braun is facilitating the management of high output stoma with B. Braun Flow Collector® which is compatible with all the Flexima® / Softima® range (High Flow products). It will reduce the need to get up at night to empty the system.

### 10 Removing the bag

Change/empty the bag if volume capacity is reached and/or if you experience odors or ballooning. Diet can affect the stool consistency, ballooning, and odor.

Pull slightly on the top of the skin protector and remove the skin protector from the top to the bottom while maintaining gentle pressure on the skin around the stoma with your other hand. Always remove the skin protector/ bag gently and carefully to reduce peristomal skin complication and additionally to reduce discomfort and pain.

Discard the ostomy bag according to the current local regulations. Do not flush down the toilet.

#### Clinical benefits & performances:

Flexima® / Softima® Active allows the collection and containment of stools or urine.

If you note any change on your stoma and/or peristomal skin, please contact your health care professional.

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