



ORTHOPAEDIC SURGERY

# AESCULAP® e.motion® Pro SYSTEM

KNEE ARTHROPLASTY  
OPERATING TECHNIQUE WITH IQ INSTRUMENTS

# AESCULAP® e.motion® Pro SYSTEM



## 1 | INTRODUCTION



The IQ e.motion® instrumentation has been designed to facilitate the workflow not only for the surgeon, but the operating room (OR) team as a whole, by enhancing ergonomics and operative efficiency. IQ stands for "Intuitive and Quick". The system offers multiple options covering different implantation philosophies that allow each surgeon to follow his/her preferred surgical technique.

- Precise and less instruments,
- quick couplings,
- ergonomic handles and
- colour coding

are some aspects that will facilitate the surgical process in the operating room.



# IQ – INTUITIVE & QUICK LESS IS MORE

The instruments as well as the instrument trays are colour coded to ease instrumentation and organization during the complete workflow:

- red = femur
- blue = tibia
- yellow = general instruments
- grey = patella

The IQ e.motion® instruments are stored in the specially developed Aesculap OrthoTray®s. Both together, the IQ instruments plus Aesculap OrthoTray® offer a high end reprocessing solution. The trays not only store the instruments in a secure and safe manner but also clearly facilitate the reprocessing procedure for the CSU (Central Sterilization Unit) as the instruments can remain in the tray during the washing process. This time saving solution generates an economic advantage and eliminates a potential source of error as complete set reassembling is needless (1).

## Aesculap Reset®

Aesculap Reset® is an intelligent improvement of the Aesculap OrthoTray® configuration. All size-specific instruments are packed such that only the sizes desired by the surgeon are used. Thus, the instrument and tray volumes in the entire instrument cycle are reduced by more than 50% (1). Aesculap Reset® facilitates, as size-specific storage and washing system, the work of all the participants in the entire process.

## NOTE

This wash tray system is only approved for the use with the cleaning validated instruments from AEscuLAP®. Complex instruments, e.g. cutting guides or instruments that are introduced in the intramedullary (IM) canal during the procedure as drills and reamers require a manual pre-cleaning according to standard requirements.

# AESCULAP® e.motion® Pro SYSTEM

## 2 | CONTENT



1	INTRODUCTION	2
2	CONTENT	4
3	INDICATIONS / PATIENT SELECTION	6
4	PREOPERATIVE PLANNING	7
5	APPROACH	8
6	ASSEMBLY INSTRUCTIONS AND INSTRUMENT HANDLING	10
7	WORKFLOW SYNOPSIS	16
8	TIBIA PREPARATION	20
	8.1 Extramedullary Alignment	
	8.2 Intramedullary Alignment	
	8.3 Tibia Resection	
	8.4 Tibia Augment Preparation	
	8.5 Tibia Wing Preparation	
	8.6 Tibia Stem Preparation	



9	<b>FEMUR PREPARATION</b>	33
	9.1 Intramedullary Femur Alignment	
	9.2 Distal Resection	
	9.3 AP Sizing and Femur Rotation	
	9.4 4-in-1 Resection	
	9.5 PS Box Preparation	
10	<b>GAP BALANCING</b>	41
	10.1 Tibia First Gap Management with Spacer	
	10.2 Tibia First Gap Management with Spreader	
	10.3 Femur First Gap Management with Spacer	
	10.4 Gap Management Strategies	
11	<b>PATELLA PREPARATION</b>	46
12	<b>TRIAL REDUCTION</b>	48
13	<b>COMPONENT IMPLANTATION</b>	50
14	<b>CEMENTING TECHNIQUE</b>	54
15	<b>INSTRUMENTS</b>	56
16	<b>SAW BLADES</b>	68
17	<b>IMPLANT DIMENSIONS</b>	69
18	<b>LOAN SYSTEMS</b>	73
19	<b>IMPLANTAT MATRIX</b>	75
20	<b>LITERATURE</b>	82



# AESCULAP® e.motion® Pro SYSTEM

## 3 | INDICATIONS/PATIENT SELECTION


### e.motion® FP

*cemented/Plasmapore®  $\mu$ -CaP*



### e.motion® UC Pro

*cemented/Plasmapore®  $\mu$ -CaP*



### e.motion® PS Pro

*cemented*

e.motion® Pro system is indicated for patients requiring primary surgery. The implant concept principle of e.motion® is based on high congruency between the femoral condyles and the mobile meniscal component and therefore requires stable collateral ligaments, medio-lateral symmetry and congruent flexion and extension gaps. The e.motion® UC Pro system is an ultra congruent PCL retaining design. The tibia does have a safety stop which allows  $\pm 30^\circ$  of rotation. The e.motion® PS Pro is a posterior stabilizing PCL sacrificing design. Due to the bone-saving design of the femoral implant, it allows implantation of a PS design without major bone loss. Another special feature is the Advanced Surface technology with its excellent gliding, wear, and allergy properties.

The successful since 2007 used surface treatment Advanced Surface, with the seven layer architecture, provides multiple advantages. The ceramic, very hard surface, with the excellent tribological properties, reduces wear by up to 60% (1, 2). At the same time, five intermediate layers are applied for the reduction of metal ion release, as an approach to prevent allergic reactions. The stability of the coating is supported by an improved elasticity and a graded degree of hardness (3, 4). For more information and contra-indications, please refer to the instructions for use TA016100.

#### NOTE

Cementless versions of e.motion implants with AS technology are available on request on custom made service.

## 4 | PREOPERATIVE PLANNING



Long leg x-ray for planning of the mechanical axis

For every Total Knee Arthroplasty, careful preoperative X-ray planning is recommended in order to determine precisely the following parameters:

- Varus/Valgus deformity
- Angle between the anatomical and mechanical femoral axes
- Entry point(s) of the intra-medullary alignment rods (manual IM technique)
- Joint line level
- Femur resection heights
- Tibia resection heights
- Component sizing
- Implant positioning
- Potential areas of bone losses and location of osteophytes

The following X-ray images are required to conduct the radiographic analysis:

- Knee joint in AP projection: knee extended, centered over the distal patella.
- Knee joint in lateral projection: knee in 30° flexion, centered above the distal patella.
- Image of the whole leg (from hip to ankle) in monopodal stance.
- Patella-tangential image (Merchant View) with the knee at 30° flexion.

The angle between the mechanical and anatomical femur axes is measured with the combination template for axis measurements. The center of the joint, the joint line and the mechanical femur axis can be measured.

To determine the tibia resection, the template showing representations of the tibial components is superimposed over and aligned with the X-ray image. The resection height is given at a 10–24 mm graduation. A complete set of radiographic templates is provided for the preoperative determination of the appropriate implant sizes. The localization of the osteophytes facilitates their removal, improving the mobility of the joint. The e.motion® knee system provides a complete set of radiographic templates in different magnitudes (PS Pro and UC Pro Tibia 1.1:1 NS416 and 1.15:1 NS417), (UC Pro Femur 1.1:1 : NE398 and 1.15:1 NE399). The results of the preoperative planning should be documented in the patient's file and available during the operative procedure for reference.

# AESCULAP® e.motion® Pro SYSTEM

5 | APPROACH




The IQ e.motion® instrumentation is designed for use with or without the use of OrthoPilot® Navigation, for both conventional and less invasive approaches to the knee.

The initial skin incision is a straight midline or slightly oblique parapatellar skin incision starting 2 to 4 cm proximal to the superior pole of the patella and extending distally to the medial aspect of the tibial tubercle. The surgeon should decide on a patient basis how long of an incision is necessary for proper visualization of the knee anatomy. A parapatellar skin incision will be of benefit to patients when attempting to kneel after the operation.


The length range of the incision is generally between 8 and 14 cm symmetrically distributed above and below the joint line. Extension of the skin incision may be necessary during the procedure depending on the patient anatomy, the soft tissues and the skin tension.

Three basic types of arthrotomies are recommended for use to carry out the intra-articular exposure: the medial parapatellar, the mid-vastus or the sub-vastus.







**Skin incision**



**Medial-parapatellar Arthrotomy**



**Mid-vastus Arthrotomy**



**Sub-vastus Arthrotomy**

### **Medial-parapatellar Arthrotomy**

With the knee in flexion or extension, the arthrotomy is performed starting proximal to the superior pole of the patella, incising the rectus femoris tendon longitudinally. Continuing the arthrotomy distally around the medial aspect of the patella, and ending medial to the tibial tubercle is then carried out.

### **Mid-vastus Arthrotomy**

With the knee in flexion, the arthrotomy is performed starting by a split of the fibers from the vastus medialis oblique (VMO), continuing distally around the medial aspect of the patella, and ending medial to the tibial tubercle.

### **Sub-vastus Arthrotomy**

With the knee in flexion, the arthrotomy is performed starting with a 4 to 6 cm incision of the fascia at the inferior border of the VMO, running horizontal to the medial aspect of the patella, continuing and ending distally medial to the medial tubercle.

### **Final exposure**

A fat pad excision is performed in order to facilitate the exposure and to improve the patella mobility. Perform the necessary medial release at this time that corresponds to the deformity. The patella can then be everted or sub-luxated laterally.




# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## 6 | ASSEMBLY INSTRUCTIONS AND INSTRUMENT HANDLING




A	TIBIA EXTRA-MEDULLARY ALIGNMENT	20
B	TIBIA INTRA-MEDULLARY ALIGNMENT	24
C	FEMUR INTRA-MEDULLARY ALIGNMENT	33
D	A/P AND ROTATION ALIGNMENT BLOC	35
E	TIBIAL-DISTAL CUTTING BLOC	34


## A | TIBIA EXTRA-MEDULLARY ALIGNMENT – ASSEMBLY INSTRUCTIONS




- Press the upper button on the bimalleolar clamp.
- Engage the support in the groove.
- When the neutral position is reached, release the button.




- Turn the wheel of the tibial alignment handle to the open position, OP-EN will be displayed.
- Engage the handle onto the bimalleolar support.
- Adjust to the neutral position.




- Push on the handle adjusting wheel to release the locking mechanism.
- Engage the holding rod in the handle.
- Release the wheel when the desired level is reached.
- Turning the wheel will allow a fine adjustment on the height.



- Engage the holding rod in one of the connection squares of the tibial cutting guide.
- Lock the assembly by turning the frontal wheel.




- The proximal fixation is set through the proximal opening of the holding rod.
- Turn the tab into a horizontal position to fix the assembly.




- The connection square of the stylus is engaged in one of the connection squares of the tibial cutting guide.
- The connection is fixed by locking the wheel on the stylus.
- The resection height is adjusted to the desired bone cut level.
- The stylus can be placed over the proximal fixation.

# AESFULAP® e.motion® Pro SYSTEM

## B | TIBIA INTRA-MEDULLARY ALIGNMENT



- Push on the button of the T-handle to release the locking-mechanism.
- Couple the T-handle to the IM rod.
- Release the button to lock the assembly.




- Choose the IM orientation sleeve corresponding to the desired posterior slope resection of the tibia (default is 0° sleeve; sleeves with 3°, 5° and 7° posterior slope are available).
- Connect the sleeve to the IM alignment system.




- Mount the assembly into the alignment block.
- Connect the alignment system to the tibia cutting guide in one of its connection squares.
- Fix the connection by locking the wheel.


## C | FEMUR INTRA-MEDULLARY ALIGNMENT



- Push on the button of the T-handle to release the locking-mechanism.
- Couple the T-handle to the IM rod.
- Release the button to lock the assembly.





- Choose the IM orientation sleeve corresponding to the desired valgus alignment (standard: 5, 6, or 7°).
- Connect the sleeve to the IM alignment system.
- Connect a distal femur contact plate (small or large).



- Mount the assembly into the alignment system.
- Connect the alignment system to the tibia cutting guide in the central connection square.
- Fix the connection by locking the wheel.


## D | A/P AND ROTATION ALIGNMENT BLOCK FOR THE FEMUR




- **Option 1:** the rotation is pre-fixed (wheel A) to a desired value before the block is put in place.
- **Option 2:** the rotation is free and the block is placed in contact with the distal femur and the posterior condyles; the rotation can be tuned by turning the posterior wheel (A), checking the alignment of the AP window (B) with the femur AP plane (Whiteside line).
- Due to the fixed distance between the pin placement holes and the anterior cortex stylus, the placed pins can be used for any femoral size chosen by the surgeon. Upsizing or downsizing the femur is achieved simply by choosing a different 4-in-1 cutting block size and placing on the same previously placed pins.

# AESCULAP® e.motion® Pro SYSTEM


## D | A/P AND ROTATION ALIGNMENT BLOCK



- The anterior point to be palpated is located on the lateral anterior cortex, avoiding the risk of anterior notching.
- If the palpation is done at the middle of the anterior femur, the grand piano sign will be bigger providing a larger surface of contact.
- The stylus can be adjusted in the caudo-cranial direction in order to get a congruence between the AP sizing and the proximo-distal sizing determined by the scale on the upper part of the stylus.



- After defining the right axial rotation of the block, if an exact femoral size is measured like in the example on the left, fix the AP sliding by tightening the corresponding screw, place 2 headless pins in the placement holes.
- By loosening the screws, and, if used, removing the posterior enhanced fixation pins, remove the orientation block.



- After defining the right axial rotation of the block, if the measured size is in between two exact sizes like in the example on the left, fix the AP sliding by tightening the corresponding screw, place 2 headless pins in the placement holes.
- By loosening the screws, and, if used, removing the posterior enhanced fixation pins, remove the orientation block.
- Choose the upper or lower size based on the assessment of the medio-lateral dimension and the flexion/extension gap situation. A smaller size will enlarge the flexion gaps; a bigger size will reduce the flexion gaps.

### NOTE

The posterior and distal thickness of the e.motion® femur differs depending between the following 3 size groups:

S = size 2, 3 = 7.0 mm; M = size 4, 5, 6 = 8.5 mm and L = 7, 8 = 10 mm.

Up- or downsizing can therefore also have an impact on the extension gap.



## E | TIBIAL-DISTAL CUTTING BLOCK

### Distal resection or tibial resection with a standard approach

- The connection to the alignment system to be used is the central one marked 'C', denoted by the green square in the left picture.
- The fixation holes for the headless pins to be used correspond to the groups marked 'C', shown by the red circles on the left picture.
- Enhanced fixation is achieved with one or two converging pins in the holes marked with the blue circles.



### Right knee tibial resection with a less invasive approach

- The connection to the alignment system to be used is the one marked "R", shown by the green square in the left picture.
- The fixation holes for the headless pins to be used correspond to the groups marked "R", shown by the red circles in the left picture.
- Enhanced fixation is achieved with one converging pin in the hole marked with the blue circle.



### Left knee tibial resection with a less invasive approach

- The connection to the alignment system to be used is the one marked "L", shown by the green square in the left picture.
- The fixation holes for the headless pins to be used correspond to the groups marked "L", shown by the red circles in the left picture.
- Enhanced fixation is achieved with one converging pin in the hole marked with the blue circle.




### NOTE



For a minimally invasive approach or when there is little space there are medialized cutting guides available as an option (see Chapter 16 Optional instruments).

# AESFULAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM




7 | WORKFLOW SYNOPSIS – TIBIA FIRST




optional




optional







6 Femur A/P sizing and rotation




7 Femur APC resections



8 Tibia keel and/or stem preparation




9 Trial reduction



10 Component implantation

OPTIONAL






PS box preparation






Patella preparation






# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

7 | WORKFLOW SYNOPSIS – FEMUR FIRST





OPTIONAL






OPTIONAL





# AESFULAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## 8 | TIBIA PREPARATION – EXTRAMEDULLARY ALIGNMENT



### 8.1 Extramedullary Referencing

- The EM alignment system assembly is placed in a parallel fashion with the frontal tibia with the leg positioned in flexion. To reach 0° tibia slope, the best way is to place two fingers between tibia and alignment system in the upper third and three fingers in the lower third.
- The bimalleolar clamp, previously set in a neutral position, is fixed around the lower limb just above the ankle joint and centered on the tibio-tarsian joint.
- Proximally, the EM alignment system can be stabilized with the proximal fixation first by engaging the longest spike between the tibia spines.

## INSTRUMENTS



Bimalleolar clamp  
NS345R



Bimalleolar clamp support  
NS344R



Alignment system handle  
NS342R




Holding rod for cutting  
guide NS341R



Tibia cutting guide  
NS334R


### Varus-valgus alignment

Pushing the knob (A) at the bimalleolar clamp, and sliding the alignment system medially or laterally allows to adjust the varus/valgus of the proximal tibia resection. The distance between the laser marked lines on the scale corresponds to a 1° adjustment for a 40 cm long tibia.



### Tibia Slope alignment

Releasing the fixation wheel (B) at the bottom part of the alignment system (by aligning OP-EN), the alignment system can be shifted anteriorly in order to increase the slope of proximal tibia resection. The distance between the laser marked lines on the scale corresponds to a 1° adjustment for a 40 cm long tibia.



---

#### NOTE

For e.motion® a slope of 0° is recommended.

---




Proximal fixation NS343R



Tibia stylus NS347R


# AESCULAP® e.motion® Pro SYSTEM

## 8 | TIBIA PREPARATION



### 8.1 Height adjustment (C)

- The resection height is determined during preoperative planning. The aim is to remove a defect on the tibial joint surface as completely as possible in order to create a level bed for the tibial plateau on intact bone.



- The planned value is set on the stylus, which is then mounted into the cutting guide. The extra-medullary alignment instrument is then lowered (wheel 3) until the stylus comes into contact with the chosen point.
- Reference to the healthy tibial plateau is helpful in determining the height of the joint line. Reference to the lowest point of the worn side of the tibia helps minimize the resection by resecting only 2 mm. Preoperative planning and surgeon preference determine which reference to use.

### NOTE

The thinnest tibial implant has a thickness of 10 mm (metal + PE) and grows in 2 mm increments.

## INSTRUMENTS



Bimalleolar clamp  
NS345R



Bimalleolar clamp support  
NS344R



Alignment system handle  
NS342R



Holding rod for cutting  
guide NS341R




Tibia cutting guide  
NS334R




Tibia stylus NS347R

- The cutting block is fixed with two headless pins in position "0". The +/-2 mm pinholes are available on the resection blocks to further adjust the resection level if needed. To avoid movements during the resection, additional pins are set in convergent holes as marked.



- The EM tibia alignment system is then disconnected from the tibia cutting guide by turning the connecting wheel counterclockwise. The proximal fixation can be removed by disengaging the spike from the tibial spine.



Headless pins 63 mm  
NP583R




Pin driver NP613R



Acculan 4 GA330 with  
drilling adapter GB664R

# AESCULAP® e.motion® Pro SYSTEM

## 8 | TIBIA PREPARATION – INTRAMEDULLARY ALIGNMENT




### 8.2 Intramedullary Alignment

- The medullary canal of the tibia is opened with the Ø 9 mm starting drill bit. The surgeon has to pay close attention of the drilling direction in order to avoid cortical violation of the posterior metaphysis.

#### NOTE

A tibial slope that is too large can result in an anterior cortical impingement if long tibial stems are used.



- After irrigation and suction, the intramedullary rod is inserted into the prepared canal using the T-handle. Once the T-handle is removed, the intramedullary alignment system is mounted on the rod with the chosen posterior slope angle sleeve (0, 3, 5, or 7°) and the cutting guide.

#### NOTE

For e.motion® a slope of 0° is recommended.

## INSTRUMENTS



Drill Ø 9 mm NE443R



T-handle NE198R



IM alignment rod NS331R



IM alignment system NS332R




Tibia cutting guide NS334R



Tibia IM stylus for orientation sleeves NS847R



- The stylus is set on the deepest point of the tibia plateau to define the 0-level cut. The height of the cut is then adjusted by turning the tuning wheel to the desired amount of resection in millimeters.






---

**NOTE**

The thinnest tibial implant has a thickness of 10 mm (metal + PE) and grows in 2 mm increments.

---

- The alignment of the cutting block can be checked with the alignment rod.




Alignment rod long  
NP471R




Tibia IM orientation  
sleeves 0°, 3°, 5°, 7°  
NS843R-NS846R

# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## 8 | TIBIA PREPARATION



- The cutting block is fixed with two headless pins in position "0". The +/-2 mm pinholes are available on the resection blocks to further adjust the resection level if needed. In order to avoid movements during the resection, additional pins are set in convergent holes.
- The cutting block is unlocked by opening the knob. The IM tibial alignment system is then removed in one step using the T-handle.



### 8.3 Tibia Resection

- Once the cutting block is positioned and fixed, the proximal tibial resection is performed. After performing the proximal tibial resection the block is removed and the resected bone taken away.
- A careful inspection of the peripheral resection is mandatory in order to check that no remaining bone stock is present. Further removal of meniscal remnants and osteophytes that encroach the posterior capsule is then performed.

#### NOTE

The protection of the surrounding soft tissues of the knee joint is paramount. A special attention has to be paid: use of Hohmann retractors, collateral retractors, PCL retractor is recommended in order to protect them during the resection.

## INSTRUMENTS



IM alignment rod  
NS331R



IM alignment  
system NS332R



Tibia IM stylus for  
orientation sleeves  
NS847R



Tibia cutting guide  
NS334R




Headless pins 63  
mm NP583R



Tibia Intramed-  
ullary orientation  
sleeve  
NS843R - NS846R



Acculan 4  
oscillating saw  
GA330






Single-use saw  
blade GE249SU  
1.27mm thick

## 8.4 Tibia Augmentation

- For bony defect on the tibia plateau, the e.motion® system offers semi tibia augments in heights of 4 mm, 8 mm and 12 mm.

After the tibia standard resection, the tibia cutting guide is removed. According to the used pin wholes on the tibia cutting block, the move block placed over the two headless pins on the side "C" or "LR" (see figure below).


Two additional headless pins are inserted in the requested depth. After removing the move block and the first two parallel pins, the tibia cutting guide can be placed over the new drilled pins.



Tibia / distal cutting guide  
NS334R



Headless pins 63 mm  
NP583R



Tibia move block  
NQ1077R




Acculan 4 GA330 with  
drilling adapter GB664R




Pin driver NP613R

# AESCULAP® e.motion® Pro SYSTEM

## 8 | TIBIA PREPARATION



When the cutting guide is fixed with two convergent pins, the hemi spacer resection can be performed. For the sagittal cut a reciprocating saw is used.



### NOTE

For tests with the trial implants the correct hemi spacer has to be clicked under the trial tibia plateau. During measuring of extension and flexion gap the height of the hemi spacer is added at the resected tibia side. An e.motion® UC/PS Pro tibia plateau which offers the possibility to screw hemi spacer has to be used for definitive implantation.

## INSTRUMENTS



Tibia/Distal cutting guide  
NS334R



Headless Pin 63 mm  
NP583R




Acculan 4  
oscillating saw GA330



Screw Driver NS410R





Tibia preparation plateau  
NS532R-NS538R



Tibia trial augment  
NS922-NS947

## 8.5 Tibia Keel Preparation

- The ideal size of the asymmetric tibia implant is determined by superposing the different tibia preparation plateau sizes onto the created surface trying to reach the best bony coverage with the proper transverse rotational alignment of the trial base plate while avoiding ML and AP overhang.
- The rotation peg is clicked into the tibia plateau. With the rotation peg the position and rotation of the meniscus component can be checked.
- The chosen tibia trial preparation is placed flush onto the tibia resection and the rotation is assessed with the help of the EM rod placed through the holder. References for the rotation are the mid-third of the anterior tuberosity and the second toe axis of the leg. These two landmarks are often not coincident with mechanical axis of the tibia and the surgeon should consider the rotation with respect to the tubercle to maintain extensor mechanism alignment. The plateau is fixed by the short headed pins in the marked holes.
- Another option is to insert the tibia and femur trial implant with the adequate trial meniscal component. The rotation peg insert of the corresponding size group helps to main the meniscal component in place. By exercising flexion extension movements combined with slight rotational stresses, the tibia plateau will find a natural position under the femur trial. This position is marked anteriorly using the electric cautery right where the plateau has a central anterior laser marking. Care should be taken to assess the stability of the extensor mechanism before accepting this 'free float' alignment of the tibial base plate.




Acculan 4 GA330  
with drilling adapter  
GB664R



Acculan 4  
oscillating saw  
GA330



Tibia preparation  
plateau NS732R-  
NS738R



Tibia trial/prep.  
plateau holder  
NQ378R



Headed pins 30 mm  
NP585R




Pin driver NP613R




rotation peg NS739R

# AESCULAP® e.motion® Pro SYSTEM

## 8 | TIBIA PREPARATION



- The holder is removed. The guiding tower of the corresponding size group is placed on the tibia plateau by engaging the posterior teeth first. The anterior part is stabilized by positioning the tower over the headed fixation pins.
- The drill with stop is first used to prepare the bone for the winglet chisel.



- The wing stem preparation is performed by using the winglet chisel corresponding to the chosen size group connected to its handle. The chisel is hammered in until flat on the tibia plateau. If no stem is used, the chisel wing stays in place as a trial. The two lever are pushed for removal of the handle. If a stem is used, the chisel is removed using the hammer.


## INSTRUMENTS



Tibia trial/preparation plateau NS532R-NS538R



Headed pins 30 mm NP585R



Guide for winglet chisel NS527R-NS529R



Drill with stop NS521R-NS523R



Acculan 4 GA330 with drilling adapter GB664R





## 8.6 Tibia shaft preparation

In case of poor bone quality, the primary fixation can be enhanced by using a stem extension. According to the surgeon's philosophy, a cemented stem or a cementless stem can be chosen.

### Option 1: priority to the tibia resection

In this case, the tibia preparation is performed following the steps described previously (§ 8.1 to § 8.4). At the last stage, instead of using the standard Ø 14 mm drill, a long drill is used for preparing the site of the future stem.

Length and diameter of this long drill should be assessed on the pre-operative X-rays. The drilling is performed through inserts for the guiding tower and the diameter (Ø 12, 14 or 16 mm) corresponds to the trial stem diameter. Three laser markings are available on the drill in order to define the right depth for short, middle or long stems. For the winglet preparation, the corresponding trial tibia stem is connected to the winglet chisel for the final preparation. Please note that this option is indicated for cemented stems.



#### NOTE

The implant stems have diameters Ø 10, 12, 14 and 16 mm in order to manage a 1 mm cement mantle thickness around the stems. For the cement mantle of the 16 mm stem the 18 mm reamer should be used.

#### NOTE

For safety reasons the long stems can only be prepared with reamers not with drills.



Obturator for tibia broach  
NS363R



Winglet chisel / trial keel  
NS791R



Osteodenser holder  
NS520R



Tibia drill sleeve for  
cemented stem  
NS547R-NS549R




Drill for cemented stem  
NS544R-NS546R



Trial stem cemented  
NE094T-NE097T,  
NE114T-NE117T,  
NE124T-NE127T

# AESCULAP® e.motion® Pro SYSTEM

## 8 | TIBIA PREPARATION




### Option 2: priority to the extension stem fixation

In this case, the medullary canal of the tibia is opened according to the preoperative planning (entry point) with the Ø 9 mm drill. The thinnest reamer is then coupled to the T-handle and inserted into the tibia medullary canal as deep as possible until a primary stability is achieved and a depth laser marking reaches the estimated level of the tibia resection (short or long stem). If not, a thicker diameter is used until stability is achieved. Once the T-handle is removed, the intra medullary alignment system is mounted on the reamer with the 0° angle sleeve (angled sleeve for slope is not possible here!) and the cutting guide. The stylus is set on the deepest point of the tibia plateau to define the 0-level cut.

### NOTE

For safety reasons, the 132 mm stems can only be prepared with the reamers and not with the drills.




The height of the cut is then adjusted by turning the tuning wheel. The alignment of the cutting block can be checked with the EM alignment rod. The cutting block is fixed with two headless pins in position "0"; the +/-2 mm pinholes are available on the resection blocks to further adjust the resection level if needed. In order to avoid movements during the resection, additional pins are set in convergent holes if necessary. The IM tibia alignment system is removed in one step with the T-handle after unlocking the cutting block from the alignment system. Please note that this option is indicated for cementless stems and the surgeon must take into account the alignment of the tibia as directed by the cementless stem since it may not coincide with the mechanical axis of the tibia.

## INSTRUMENTS



Reamer for cementless stem NE154R-NE158R



IM alignment rod NS331R




IM alignment system NS332R



Tibia IM stylus for orientation sleeves NS847R



Tibia cutting guide NS334R




0° orientation sleeve NS843R

## 9 | FEMUR PREPARATION

### 9.1 Femur Intra-medullary Alignment

- The medullary canal of the femur is opened according to the preoperative planning (entry point) with the drill  $\varnothing$  9 mm. The rod is inserted into the intra medullary canal using the T-handle. Once the rod is inserted, the T-handle can be removed.
  - For evaluation of the appropriate distal cut, the femur size group is measured. The final size can be decided later.
  - The femur sizing is achieved by reading frontally the marked size on the scale when the stylus tip is placed at the intended exit point of the saw blade on the anterior lateral cortex in order to avoid any notching.
- 
- In order to compensate the anatomical valgus angulation of the femoral bone, the appropriate angle sleeve  $5^\circ$ ,  $6^\circ$  or  $7^\circ$  according to the preoperative planning is set into the intra-medullary alignment system (angle sleeve  $8^\circ$  and  $9^\circ$  are available on demand). The distal femur contact plate and the cutting block are connected to this system. The assembly is placed on the IM rod in contact with at least one distal condyle.
  - The planned height of the distal resection is adjusted by turning the wheel (A) until the desired thickness matches the anterior laser marking. The standard resection corresponds to the distal thickness of the implant and is 7 (size group S), 9 (size group M) and 10 mm (size group L) depending on the size group.



Femur alignment block NS580R



Drill  $\varnothing$  9 mm NE443R



Acculan 4 GA330 with drilling adapter GB664R




T-handle NE198R



Tibia alignment system NS332R



Distal femur contact plate NS333R, NS834R




Femur orient. sleeve NS335R-NS337R



Tibia cutt. guide NS334R


# AESCULAP® e.motion® Pro SYSTEM

## 9 | FEMUR PREPARATION



### 9.2 Distal Resection

- The cutting block is fixed with two headless pins in position "0". To avoid movement during resection, additional pins are set in convergent holes.




- The intra-medullary alignment system is completely removed in one step with the T-handle by unlocking the connection to the cutting guide.
- The distal femoral resection is performed by sawing through the slot with a 1.27 mm thick oscillating saw blade. Make sure that the resection is fully completed and that no remaining bone structures are prominent to the resection plane.
- Pins and cutting block are removed.

### NOTE

Please always pay a great care to the lateral structures by protecting them if necessary by the use of Hohmann retractors.

## INSTRUMENTS



IM alignment rod  
NS331R



Tibia alignment  
system NS332R




Distal femur contact  
plate NS333R,  
NS834R



Femur orient. sleeve  
NS335R-NS337R



Tibia cutt. guide  
NS334R




Headl. pins 63 mm  
NP583R



Acculan 4 GA330  
with drilling adapter  
GB664R

### 9.3 Final Femur Sizing and Rotation

- The ML size of the resected femur should be checked with the ML femoral sizing gauge.
- The femur alignment block is placed flush onto the resected distal surface of the femur. The posterior foot plate must be in contact with the posterior condyles. The femoral alignment block is fixed with two headless pins against the distal femur through the posterior holes.



- The femur sizing is achieved by reading frontally the marked size on the scale when the stylus tip is placed at the intended exit point of the saw blade on the anterior lateral cortex in order to avoid any notching. A scale on the surface of the stylus indicates the femur size depth and the position can then be fixed by tightening the screw.



#### NOTE

To avoid notching at the anterior femur, make sure that the stylus fixation screw on the femur alignment block is not too loose.




T-handle NE198R



Tibia protection plate  
NQ377R



Acculan 4 GA330 with  
drilling adapter GB664R



ML femoral size gauge  
NS581R




Pin driver NP613R




Headless pins 63 mm  
NP583R

# AESCULAP® e.motion® Pro SYSTEM

## 9 | FEMUR PREPARATION



- It is possible to adjust the external rotation by moving the posterior lever arm in the right direction (clockwise for right knees, anticlockwise for left knees). The rotational position is confirmed by assessing the transepicondylar axis perpendicularity or by checking the Whiteside's line through the slot at the middle of the instrument. Size and rotation are fixed by tightening the screw at the bottom lever arm.
- The insertion of standard fixation pins on the medial and lateral aspect of the femur alignment block facilitates referencing of the epicondyles.



- Two long headless pins are fixed through the 2 frontal holes in order to reference the position of the 4-in-1 cutting guide. It is recommended to check the level of the anterior resection by using the check plate in the alignment block slots. The size to choose is to be read on the scale (see § 6 handling instructions).
- The posterior pins and the block are removed, leaving the headless pins in place.

### INSTRUMENTS



Femur alignment block  
NS580R



Headless pins 63 mm  
NP583R



Pin driver NP613R



Acculan 4 GA330 with  
drilling adapter GB664R




#### 9.4 Femur Anterior, Posterior and Chamfer Resections

- The 4-in-1 cutting guide that matches the femur size is placed over the two headless pins into the marked "0" mm pinhole and pressed onto the distal resection. It is advised to check the level of the anterior resection by using the check plate in the alignment block slots before placing the converging headed pins for fixation.
- Before fixing the guide with convergent headless pins, it is possible to adjust the AP position by using the holes marked +/- 1.5 mm in order to remain as close as possible to the anterior cortex without notching it.



- The resections are performed as follows: anterior cut, posterior cut, removal of sizing pins, posterior chamfer, anterior chamfer. Thereby, the maximum distal contact surface and cutting block fixation is preserved up to the last resection, ensuring stability.
- Convergent pins and cutting guide are removed, and the resections are carefully checked in order to detect any remaining bone stock.



4-in-1 femur cutting guide NS582R-NS588R




Cutting depth check blade NS850R




Acculan 4 oscillating saw GA330

# AESCULAP® e.motion® Pro SYSTEM

## 9 | FEMUR PREPARATION



- The quality of the resections and the fit of the prosthesis can be assessed by placing the femur trial implant onto the bone preparation. Using the femur holder, make sure to apply a force toward anterior in order to avoid a flexed position.
- For reposition of the joint, the trial meniscus and the trial tibia plateau are inserted.



- For downsizing the femur, a smaller 4-in-1 cutting guide is placed directly onto the same anterior headless pins using the same holes as previously (-1.5/0/+1.5). Since the reference is anterior, you will achieve the same anterior cut but recut the posterior condyles, the posterior chamfer as well as the anterior chamfer. This will open the flexion gap correspondingly.

## INSTRUMENTS



Insert for NS600R,  
NS601-NS603



Holder for trial  
femurs NS600R



Trial femurs NS740K-  
NS759K (PS Pro),  
NE702K-NE708K,  
NE752K-NE758K,  
NS981K-NS986K (UC)



Tibia trial plateau  
NS732R-NS738R



Tibia trial holder  
NQ378R




Trial meniscus  
NS772-NS778,  
NS782-NS788 (PS Pro),  
NS622-NS658 (UC)




Complementary  
plates NS675-NS697

## 9.5 PS Box Preparation

- The femur pegs are drilled before box preparation.
- The trial femur implant is placed onto the prepared femur using the corresponding holder making sure to apply an extension force anteriorly in order to avoid a flexed position. The trial femur implant is fixed along the proximal trochlear groove with two headed pins.



- The miller guide is fixed into the holes of the distal trial femur. An additional pin avoids a movement of the miller guide.



Femur box milling guide NS797R-NS799R



Femur drill peg NE458R



Headed pin 50 mm NP586R



Pin driver NP613R



Acculan 4 GA330 with drilling adapter GB664R




Headless pins 63 mm NP583R



Drill 9 mm NE443R

# AESCULAP® e.motion® Pro SYSTEM


## 9 | FEMUR PREPARATION



- Use the special box miller for the chosen size group to prepare the small box.
- Be careful with the soft parts of the posterior femur.
- Use a protection plate for the tibia.

### NOTE

The area between posterior capsule and intercondylar notch should be coagulated in order to prevent postoperative uncontrolled bleeding.



- The rest of the bone at the medial and lateral inner box wall is removed with a saw blade or a flat chisel. In case of chisel use the two slots on the miller guide for the chisel. Care must be taken not to cut too deep into the femur.
- The miller guide is removed and the trial femur implant can be used with trial PS Pro meniscal components.
- Optionally the trial can be inserted for a simulation of the camp-post mechanism.

## INSTRUMENTS



Trial femurs NS740K-NS759K (PS Pro), NE702K-NE708K, NE752K-NE758K, NS981K-NS986K (UC)



Femur box milling guide NS797R-NS799R



Femur box miller NS794R-NS796R



Trial cam NS763R-NS767R



Acculan 4 oscillating saw GA330




Handle for trial cam NS761R

**10.1 Tibia First – Measurement with Spacers**


After performing the tibia resection, check the plane of the resection by inserting the thinnest spacer block (10 mm) in the joint. If the resection needs correction then apply the cutting block accordingly and recut the proximal tibia. The soft tissue gaps can be assessed by applying a varus/valgus stress in extension and in flexion. If the joint is too lax, insert the next spacer and repeat the operation until a spacer thickness allows the knee to reach a stable point in flexion and extension.

**NOTE**


The PCL must be released and removed prior to assessing gaps in flexion and extension since it will increase the flexion gaps once removed.



- If the medial and lateral gaps are asymmetrical, it is necessary to perform the appropriate release on the contracted side and then repeat the gap measurements with the spacers until stability is reached.
- If the flexion and extension gaps are incongruent then please refer to the chapter 10.4 strategies and define the right corrective action.
- The thickness of the last spacer that allows good balance and stability of the knee corresponds to the needed polyethylene thickness that should be used.
- At each step, the leg axis can be checked by inserting the alignment rod through the spacer handle; the rod should point respectively at the femoral head center and the ankle joint center.
- The measurements can also be done after the distal resection is performed by adding the distal cut spacer of the corresponding size group (S, M or L) for the extension measurement.



Spacer NS852R-NS854R




Control rod long NP471R



Complementary Spacer NS497 - NS499

# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

10 | GAP BALANCING



## 10.2 Optional Tibia First – Measurement with Distractor


- After performing the tibia resection, check the plane of the resection so that it correspond to the mechanical axis of the tibia. Insert the distractor into the joint and use the clamp to distract sequentially the medial and lateral gaps in extension.
- If the medial and lateral gaps are asymmetrical, it is necessary to perform an appropriate release on the contracted side and then repeat the gap measurements.

---

### NOTE

The distractor set is included in the navigation instrument set (NP138).

---



- When the joint is balanced in extension, note the thickness of the gaps, and move to the flexion gap measurement and repeat the same operation. In flexion, the possible future rotation of the femoral component should be taken into account.
- When the flexion gaps (FG) differ from the extension gaps (EG), calculate the needed thickness of the distal resection in order to equalize flexion and extension:  
distal resection = 9 mm - EG + FG.

---

### NOTE

The PCL must be released and removed prior to this step since its removal will increase the flexion gaps.

---

## INSTRUMENTS



Distraction forceps  
NP609R



Femur-tibia distractor  
NP604R




### 10.3 Femur First – Measurement with Spacers


- After completion of the femoral and tibial resections, the trial femur implant is placed on the femur. The height of the resection and flexion/extension gaps can be checked by inserting the spacers like in chapter 10.1.




- Alternatively, to measure without femur trial, the height of the distal cut spacer is clicked on the gap spacer. Now the gaps in flexion and extension can be checked.



Tibia cut spacer  
NS852R-NS854R



Added femur cut spacer  
NS497-NS499






Alignment rod long  
NP471R

# AESCULAP® e.motion® Pro SYSTEM







10 | GAP BALANCING

## 10.4 Strategies

When the flexion and extension gaps are incongruent, an individualized strategy has to be defined in order to correct it. The table presents some possible options to follow in order to correct a situation where the flexion and extension gaps are not both equally optimal but either tight or wide. This does not claim to be exhaustive of all options for gap management.

		optimal
Extension gap	optimal	
	small	 <ul style="list-style-type: none"> <li>▪ posterior capsule release</li> <li>▪ increase distal cut</li> </ul>
	wide	 <ul style="list-style-type: none"> <li>▪ decrease distal cut</li> <li>▪ downsize femur and thicker insert</li> </ul>

## Flexion gap

	small	wide
	<ul style="list-style-type: none"> <li>▪ increase tibia slope</li> <li>▪ downsize the femur</li> </ul>	
	<ul style="list-style-type: none"> <li>▪ thinner insert</li> <li>▪ increase tibia cut</li> </ul>	
	<ul style="list-style-type: none"> <li>▪ downsize femur and thicker insert</li> <li>▪ downsize femur and decrease distal cut</li> <li>▪ decrease distal cut</li> </ul>	

### Overview femur/tibia compatibility

Size	F2	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
T2	Standard combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination
T3	Standard combination	Standard combination	Possible combination	Possible combination	Standard combination	Standard combination	Possible combination	Possible combination	Possible combination	Possible combination
T4	Possible combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Possible combination	Possible combination
T5	Possible combination	Possible combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination
T6	Possible combination	Possible combination	Possible combination	Possible combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination	Standard combination
T7	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Standard combination	Standard combination	Standard combination	Standard combination
T8	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Possible combination	Standard combination	Standard combination

Tibia > Femur = no limitation

Femur > Tibia = maximal two femur sizes larger are allowed

■ Standard combination

□ Possible combination


**NOTE:** The meniscal component size is based on femoral size. The same combinations apply for the Narrow (N) components. All patella sizes can be combined with any femur size.

# AESCULAP® e.motion® Pro SYSTEM


## 11 | PATELLA PREPARATION



- The thickness of the patella is measured using the caliper. This thickness should not be exceeded after implantation of the patella implant. The level of bone resection is calculated. A minimum thickness of remaining patella bone should be no less than 12 mm.
- The patella is clamped and the level of the resection is adjusted by turning the resection depth wheel (A) to the planned level of remaining patellar bone thickness.
- The resection is performed through the cutting slot with a 1.27 mm thick saw blade.



## INSTRUMENTS



Caliper AA847R





Patella resection clamp  
NS840R



Acculan 4  
oscillating saw GA330

- The patella resection clamp is removed. The patella drill/impaction clamp is set onto the osteotomized patellar surface choosing a medialized position to recreate the resected apex of the articular surface; the trial patella can be placed on top of the drill guide in order to check its position to the medial rim and appropriate positioning in the superior and inferior direction.
- The pegs of the implant are drilled through the holes with the  $\varnothing$  6 mm drill until the stop is reached. The size of the patella is established with the corresponding trial patellar implant.




Patella drill/impaction clamp NS841R



Acculan 4 GA330 with drilling adapter GB664R




Drill with stop  $\varnothing$  6 mm NQ449R




Trial patella NQ281-NQ285

# AESCULAP® e.motion® Pro SYSTEM

12 | TRIAL REDUCTION



- The trial femoral and tibial implants are placed onto the prepared bony surfaces.
- The meniscus trial corresponding to the gap measurements with the spacer or the distractor is placed between both trial implants.
- The same meniscus trials are used for left and right side. The main upper part of the trial design corresponds to the final design together with the complementary plates the desired height is achieved. With the meniscus trial the medialized center of rotation can be simulated. The R and L on the underside of the complementary plate indicates for which side of the joint the connection is suitable.
- The stability of the joint is assessed by applying varus/valgus stresses in extension and flexion. If the joint appears to be lax (opening of gaps under stress), then a thicker trial meniscus surface is tested.
- The range of motion is assessed. Intra-operative limited extension and flexion and marked hyperextension must be avoided.



## TIP

If the PS Pro version is used, check in hyperextension if the post does not impinge with anterior bone. Remove bone if necessary.


## NOTE

In case tibia augments have been prepared, they need to be added to the tibia plateau.


## INSTRUMENTS



Tibia trial/preparation plateau NS732R-NS738R



Tibia trial/prep. plateau holder NQ378R



Trial gliding surface NS772-NS788 (PS Pro), NS622-NS650 (UC)



Trial spacer 4, 6, 8 mm NS675-NS697



Trial rotation peg NS739R





# AESFULAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## 13 | PREPARATION AND ASSEMBLY OF EXTENSION STEMS



- For the assembly of the extension stem on the final implant the stem has to be tightened with a torque of 20 Nm. It is recommended to tighten the extension stem on the table and ensure that the components are hold by an assistant during the tightening.



- The rotation peg is tightened with the torch wrench for 10 Nm.

## INSTRUMENTS



Center torque  
for stem fixation  
NS570R



Torque wrench  
NE184RM



Stem adapter for  
NE184RM  
for extension stems  
Ø 12, 14 mm NE185R



Stem adapter for  
NE184RM  
for extension stems  
Ø 10 mm NS835R


- If the peek plug is used as an alternative to the opturator to close the tibia plateau, use the screw driver to screw it into the plateau.

The following implant sequence is recommended:

- Tibia implant
- Femur implant
- Meniscal component
- Patella implant

#### NOTE


If an opturator is used, it need to be tightened with the key (NS378R).



#### OPTION

The rotation peg can also be assembled to the tibia implant after the cement has hardened.

- The final tibia implant is brought precisely into the predefined position. The final positioning is achieved with the help of the tibia impactor. The meniscal component can be assembled to the tibia implant before implantation. With the torque wrench plus adapter and a counter holder the 10 Nm can be applied to the assembly.



Tibia plateau impactor NS425



Tibia implant and rotation peg



Torque wrench 10 Nm NE160R




Adapter for torque wrench NP450R  
Ø 4.5 mm



Peek plug NN260P




Screw driver SW 3.5 NS423R




Stem tightening key NS378R

# AESCULAP® e.motion® Pro SYSTEM

## 13 | COMPONENT IMPLANTATION



- Using the femur holder and the insert of the corresponding size group, the final femur implant is brought into alignment and implanted. Care must be taken to assure the holder is properly seated and attached to the femoral implant so that it does not dislodge during cementing. A special attention has to be paid to the sagittal orientation: forcing the holder to the anterior direction helps to avoid an implantation in a flexion position.
- The femur holder is opened by turning the handle counter-clockwise.




- The femoral impactor is used to knock the implant into place.


## INSTRUMENTS



Implant holding/insertion instrument NS600R



Femur insert to NS600R, NS601-NS603




Femur impactor NS424



Femur implant

- The meniscal component is placed over the rotation peg (UC or PS Pro).



**NOTE**

It may be prudent to use a trial insert and recheck joint motion and stability after the cement has cured before deciding on the final type and thickness of the meniscus insert. Therefore the trial rotation peg (NS540P) is screwed into the final tibia implant.

The patella is implanted using the patella impaction clamp and the concave plastic cap, which allows good transmission of forces during the cement hardening process.



**NOTE**

PS prostheses have a longer trochlea cutout for technical reasons. The patella tracking must therefore be checked. In the case of patella clunk syndrome, the patella must be modified or replaced with an artificial patella.



Gliding surface




Patella drill/impaction clamp NS841R



Inlay for NS841R, NS842



Patella implant



Trial rotation peg NS540P

# AESCULAP® e.motion® Pro SYSTEM

## 14 | CEMENTING TECHNIQUE

- Regardless of what fixation method is utilized it is critical that correct techniques are employed in order to avoid complications and early failure. Also, even with accurate cuts it is important to ensure that components are fully seated, as it is easy for this to be obscured when cementing is taking place. Varus-valgus alignment can be significantly affected by unequal medial-lateral cement mantles and poorly seated components and there can be a tendency to place femoral components in relatively flexed positions if specific care is not taken.
- It should also be noted that when definitive components are cemented in, they may prove more stable and seat better than the trials, which are often a little loose. It is therefore worthwhile to recheck the balancing and stability at this point so that further adjustments can be made if necessary. It has been possible to relate poor cementing techniques to early and continuous component migration, which in turn is of positive prognostic significance when predicting aseptic loosening so proper attention to the cementation steps must be taken.
- Preparation of the bony surfaces and cancellous bone should be performed with pulsatile type lavage with the knee under a pressure tourniquet. This step allows for optimal cement penetration and interlocking to the bony prepared surfaces and also removes bone debris that can serve as third body particles that increase polyethylene wear after surgery. The surfaces should be properly dried prior to cementation and appropriate exposure of all bony surfaces achieved. All of the surfaces should be pressurized for optimal cement penetration. Emphasizing the importance of effective cementation of the posterior femoral condylar surfaces is also recommended since it can have a significant effect on the longevity of the fixation of the femoral implant. A further point worth noting is that if holding the knee out in full extension while cement is hardening is used to compress components down and possibly improve cement intrusion.
- Care should be taken to completely remove all excess cement that protrudes from the implant bone interface. Any remnants of overhanging cement can impinge on surrounding soft tissue or can provide a source of debris that can serve as a generator of third body wear and may contribute to the demise of the fixation earlier than expected.







- After cement polymerization and removal of all cement excess, thoroughly irrigate the joint. If a tourniquet is used, hemostasis is achieved after its deflation. Close soft tissue in the normal layered fashion.

# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

15 | INSTRUMENTS



## SETS

Item No.	Description	Container recommended	Lid	Height of tray incl. lid
<b>NS760</b>	<b>IQ e.motion<sup>®</sup> Pro Instrumentation</b>			
<b>Consisting of:</b>				
NS701	IQ e.motion <sup>®</sup> Set General Instruments	JK444	JK489	119 mm
NS702	IQ e.motion <sup>®</sup> Set Tibia-Instruments	JK444	JK489	119 mm
NS703	IQ e.motion <sup>®</sup> Set Femur Preparation	JK442	JK489	89 mm
NS856	IQ e.motion <sup>®</sup> Pro Set Tibia Preparation	JK444	JK489	119 mm
	Insert + UC Trials for 856			
NS706	IQ e.motion <sup>®</sup> UC Set Femur Trial Implants	JK444	JK489	119 mm

Item No.	Description	Container recommended	Lid	Height of tray incl. lid
----------	-------------	-----------------------	-----	--------------------------

**Additional sets needed for e.motion® PS Pro**

**Consisting of:**

NS770	IQ e.motion® PS Pro Set Femur Trial Implants + Preparation	JK444	JK489	119 mm
NS780	IQ e.motion® PS Pro Set Trial Meniscal Components	JK441	JK489	69 mm

**OPTIONAL SETS**

Item No.	Description	Container recommended	Lid	Height of tray incl. lid
----------	-------------	-----------------------	-----	--------------------------

NS858	IQ e.motion® Set Tibia Extension Set	JK442	JK489	89 mm
NS768	IQ e.motion® PS Pro Set Femur Trial Implants standard sizes	JK444	JK489	119 mm

**e.motion® Pro Standard/Pro Tibia Augments**

NS910	IQ e.motion® Set Tibia-Augments	JK342	JK389	89 mm
-------	---------------------------------	-------	-------	-------

**e.motion® Pro Femur Trial Implant, Set is identical to NS706, but for F4-F6 in "Narrow" version.**

NS908	IQ e.motion® FP/UC Set Femur trials incl. narrow sizes (like NS706, but trials 4, 5, 6 in narrow version)	JK444	JK489	119 mm
-------	---	-------	-------	--------

**Patella**


NS709	IQ Set Patella Instruments	JK444	JK489	119 mm
-------	----------------------------	-------	-------	--------

**Navigation**


NP138	IQ Set Navigation Instruments	JK444	JK489	119 mm
-------	-------------------------------	-------	-------	--------

# AESCULAP® e.motion® Pro SYSTEM

## NS701 | GENERAL INSTRUMENTS




## NS702 | TIBIA-INSTRUMENTS




The lower part of the tray is empty and can be used for optional instruments. In case the mini stem (NB100K) is used, the torque wrench (NE184RM) with adapter (NE185R) can be stored here. It is not necessary to open the extension stem tray NS858.

# NS856 | TIBIA PREPARATION




# UC TRIAL MENISCAL COMPONENTS






# AESCU LAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## NS703 | FEMUR PREPARATION




## NS706 | FEMUR TRIAL IMPLANTS W/O NARROW FEMUR TRIALS






## NS770 | PS PRO FEMUR PREPARATION (F4-F6 IN NARROW)




## NS780 | PS PRO MENISCAL COMPONENTS




# AESFULAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## NS858 | TIBIA-EXTENSION STEMS




## NS768 | PS PRO TRIAL FEMUR STANDARD SIZES




**TF075**  
Graphic template for NS769R  
**JA415R**  
Lid for tray

## NS910 | TIBIA-AUGMENTS




## NS908 | FEMUR TRIAL IMPLANTS W/ NARROW FEMUR TRIALS



# AESCULAP® e.motion® Pro SYSTEM

NS709 | PATELLA PREPARATION



# NOTES

A series of horizontal dotted lines for taking notes, spaced evenly down the page.



# AESCULAP® e.motion® Pro SYSTEM

15 | OPTIONAL INSTRUMENTS


## GENERAL



NP609R Gap distractor for NP604R




NP604R Femur-tibia gap measuring gauge



NM640 Force controlled spreader set



NE150R Leg positioner for TKA  
NE153R Fixation Frame



Pin set (NP742R, NP743R,  
NP748R, NP749R, NP750R)

## FEMUR



NS578R Femur orientation sleeve 8°



NS579R Femur orientation sleeve 9°



## TIBIA



NS406R medialised cutting guide left



NS407R medialised cutting guide right



NS861R FGT tibia correction bloc 2° var/val



NS879R IQ FGT counter guide for NS861R



NS863R FGT tibia EM alignment system (Assembly instructions in O47302)



NE425R tibia stylus




NS844R Tibia-IM-orientation sleeve 3°



NS845R Tibia IM orientation sleeve 5°




NS846R Tibia IM orientation sleeve 7°




NP684R Slap hammer

## STORAGE OPTIONAL INSTRUMENTS



NQ1429R tray optional instruments large, lid JA455R



NE1029R tray optional instruments small, lid JA415R

### NOTE









For the optional trays the following containers and lids are recommended:

NQ1429R: Container JK442, Lid JK489



NE1029R: Container JK342, Lid JK389

# AESCULAP® e.motion® Pro SYSTEM

16 | SAW BLADES

System	Item No.	Width	Thickness	Saw blades 
AESCULAP® Acculan 3 Ti, Acculan 4 Length 75 mm	GE231SU	9 mm	1.27 mm	
AESCULAP® Acculan 3 Ti, Acculan 4 Length 90 mm	GE233SU	13 mm	1.27 mm	
	GE236SU	13 mm	1.27 mm	
	GE241SU	19 mm	1.27 mm	
	GE246SU	23 mm	1.27 mm	
AESCULAP® Acculan 3 Ti, Acculan 4 Length 100 mm	GE249SU	19 mm	1.27 mm	
Stryker System 4-7 Length 90 mm	GE330SU	13 mm	1.27 mm	
	GE331SU	19 mm	1.27 mm	
	GE332SU	25 mm	1.27 mm	
Synthes Trauma Recon System Battery Power Line Battery Power Line II Length 90 mm	GE323SU	13 mm	1.27 mm	
	GE326SU	25 mm	1.27 mm	
Zimmer-Biomet Universal Length 90 mm	GE327SU	13 mm	1.27 mm	
	GE329SU	25 mm	1.27 mm	

A complete overview of all saw blades with Aesculap® couplings are listed in our Burrs & Blades catalog: 017599.

System	Saw blades for reciprocating saws 75/10/1.0/1.2 mm	Saw blade for reciprocating saws 75/12/1.0/1.2 mm
Acculan 3 Ti, Acculan 4	 GC769R	 GC771R


## 17 | IMPLANT DIMENSIONS

### UC PRO FEMUR

AP-/ML-Dimensions [mm] of the e.motion® UC Pro/FP femoral implants for necessary application of intra medullary nail

Measurements in mm


Size	AP	ML
F2	20	18
F3	22	19
F4	24	20
F5	27	21
F6	29	22
F7	31	23
F8	33	25



### FEMUR COMPONENTS

Measurements in mm

Size	ML 2	AP	Box	Th
F2	56	50	37	7
F3	60	54	40	7
F4	64	58	43	8.5
F4 N	60	58	43	8.5
F5	68	62	46	8.5
F5 N	64	62	46	8.5
F6	72	66	49	8.5
F6 N	68	66	49	8.5
F7	76	70	52	10
F8	80	74	55	10



# AESFULAP<sup>®</sup> e.motion<sup>®</sup> Pro SYSTEM

## 17 | IMPLANT DIMENSIONS


### PS PRO FEMUR

#### Femur Components

Most important dimensions [mm] of the e.motion<sup>®</sup> PS Pro femoral implants.

Measurements in mm



Size	ML 2	AP	Box	Th
F2	56	50	37	7
F3	60	54	40	7
F4N	60	58	43	8.5
F4	64	58	43	8.5
F5N	64	62	46	8.5
F5	68	62	46	8.5
F6N	68	65	49	8.5
F6	72	66	49	8.5
F7	76	70	52	10
F8	80	74	55	10



AP-/ML-Dimensions [mm] of the e.motion<sup>®</sup> PS Pro femoral implants for necessary application of intra medullary nails.

Measurements in mm

Size	AP	ML 1
F2	30	18
F3	33	19
F4N	36	20
F4	36	20
F5N	38	21
F5	38	21
F6N	41	22
F6	41	22
F7	42	23
F8	45	24



## 17 | IMPLANT DIMENSIONS




### UC/PS TIBIA

#### Tibia Components

Most important dimensions [mm] of the e.motion® UC Pro tibia implants.

Measurements in mm

Size	ML	AP	AP/ML
T2	63	41	0.7
T3	67	44	0.7
T4	71	47	0.7
T5	75	50	0.7
T6	79	53	0.7
T7	83	56	0.7
T8	87	59	0.7




#### Patella Components

Most important dimensions [mm] of the patella implants.

Measurements in mm

Size	D	Th
1	26	7
2	29	8
3	32	9
4	35	10
5	38	11



# AESCULAP® e.motion® Pro SYSTEM

## 17 | IMPLANT DIMENSIONS

### EXTENSION STEM LENGTHS


Measurements in mm

	T0-T5
Tibia keel length	28
Tibia keel + Obturator	35
Tibia keel + Stem 12 mm	40
Tibia keel + Stem 52 mm	80
Tibia keel + Stem 92 mm	120
Tibia keel + Stem 132 mm	160

The overall length of the tibia plateau with the respective extension stem is given by the tibia keel length 28 mm and the obturator 7 mm or the stem length 12 mm, 52 mm, 92 mm or 132 mm.

### NOTE

Obturator screws and 12 mm extension stems are available in 14 mm diameter. All other extension stems are available in diameters 10, 12 and 14 mm.





## IMPLANTS

Item No.	Description
LSET-K0075	e.motion® UC/FP Femur R
LSET-K0076	e.motion® UC/FP Femur L
LSET-K0021	e.motion® UC/FP Femur R cementless
LSET-K0026	e.motion® UC/FP Femur L cementless
LSET-K0150	e.motion® PS Pro Femur L
LSET-K0151	e.motion® PS Pro Femur R
LSET-K0161	e.motion® UC Pro Tibia L
LSET-K0162	e.motion® UC Pro Tibia R
LSET-K0155	e.motion® UC Pro Meniscuscomp. L
LSET-K0156	e.motion® UC Pro Meniscuscomp. R
LSET-K0163	e.motion® PS Pro Meniscuscomp. L
LSET-K0164	e.motion® PS Pro Meniscuscomp. R
LSET-K0131	AESULAP® Tibia stems
LSET-K0210	e.motion® Pro Tibia Augments
LSET-K0176	AS e.motion® PS Pro Femur + Tibia L
LSET-K0177	AS e.motion® PS Pro Femur + Tibia R
LSET-K0178	AS e.motion® UC Pro Femur + Tibia L
LSET-K0179	AS e.motion® UC Pro Femur + Tibia R
LSET-K0211	AS e.motion® Pro Tib. Aug.
LSET-K0132	AS AESULAP® Tibia Stems
LSET-K0041	e.motion® Patella

# AESCULAP® e.motion® Pro SYSTEM

18 | LOAN SYSTEMS

## INSTRUMENTS

Item No.	Description
LSET-K0154	IQ e.motion® UC Pro Basic
LSET-K0148	IQ e.motion® PS Pro Extension to K0154
LSET-K0198	IQ e.motion® FGT
LSET-K0149	IQ e.motion® Pro Tibia Stem preparation
LSET-K0209	IQ e.motion® Tibia Augment preparation

## OPTIONAL INSTRUMENTS

Item No.	Description
LSET-K0165	IQ e.motion® UC Pro Tibia (only tibia preparation)
LSET-K0051	IQ Navigation
LSET-K0130	IQ Patella

## Aesculap Reset® INSTRUMENTS

Item No.	Description
LSET-K0199LT	Aesculap Reset® IQ e.motion® PS Pro navigated
LSET-K0206	Aesculap Reset® IQ e.motion® UC Pro
LSET-NS500LT	Aesculap Reset® IQ e.motion® Pro Basic
LSET-NS501LT	Aesculap Reset® IQ e.motion® Pro manual
LSET-NS509LT	Aesculap Reset® IQ e.motion® Pro Tibia-Prep.
LSET-NS912LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F2
LSET-NS913LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F3
LSET-NS914LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F4
LSET-NS915LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F5
LSET-NS916LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F6
LSET-NS917LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F7
LSET-NS918LT	Aesculap Reset® IQ e.motion® PS Pro Fem-Prep. F8

## X-RAY TEMPLATES

Item No.	Description
NS416	PS Pro Femur – PS Pro/UC Pro Tibia – Standard Stems 1.10:1
NS417	PS Pro Femur – PS Pro/UC Pro Tibia – Standard Stems 1.15:1
NE398	FP/UC Femur, Patella 1.10:1
NE399	FP/UC Femur Patella 1.15:1

## e.motion® Pro Implant Matrix – Femoral parts

**Patella****Types:** F2–F8

P1 NX041

P2 NX042

P3 NX043

P4 NX044

P5 NX045

**Femur FP/UC cemented**

Types:	F2	F3	F4	F4N	F5	F5N	F6	F6N	F7	F8
Left	N0502K	N0503K	N0504K	N0817K	N0505K	N0818K	N0506K	N0819K	N0507K	N0508K
Right	N0602K	N0603K	N0604K	N0917K	N0605K	N0918K	N0606K	N0919K	N0607K	N0608K

**Femur FP/UC cementless**

Types:	F2	F3	F4	F4N	F5	F5N	F6	F6N	F7	F8
Left	N0582K	N0583K	N0584K	N0837K	N0585K	N0838K	N0586K	N0839K	N0587K	N0588K
Right	N0682K	N0683K	N0684K	N0937K	N0685K	N0938K	N0686K	N0939K	N0687K	N0688K

**Femur PS Pro cemented**

Types:	F2	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NX700K	NX701K	NX702K	NX703K	NX704K	NX705K	NX706K	NX707K	NX708K	NX709K
Right	NX750K	NX751K	NX752K	NX753K	NX754K	NX755K	NX756K	NX757K	NX758K	NX759K

# AESCULAP® e.motion® Pro SYSTEM

19 | IMPLANT MATRIX

## e.motion® Pro Implant Matrix – Tibial implants



### Tibia UC Pro/PS Pro cemented

Types:	T2	T3	T4	T5	T6	T7	T8
Links	NX732K	NX733K	NX734K	NX735K	NX736K	NX737K	NX738K
Rechts	NX782K	NX783K	NX784K	NX785K	NX786K	NX787K	NX788K



### Tibia-Obturator for UC Pro/PS Pro

Types:	Ø 14 mm
Standard	NN264K
AS	NN264Z



### Tibia short stem cementless

Types:	Ø 14 mm
Länge mm	12
Standard	NB100K
AS	NB100Z



### PEEK Plug

Ø 14 mm
NN260P

## e.motion® Pro Implant Matrix – Tibial implants



### e.motion® UC Pro/PS Pro Tibia-Augments Medial

4 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left + Right	NX602K	NX603K	NX604K	NX605K	NX606K	NX607K	NX608K

8 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX622K	NX623K	NX624K	NX625K	NX626K	NX627K	NX628K
Right	NX632K	NX633K	NX634K	NX635K	NX636K	NX637K	NX638K

12 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX662K	NX663K	NX664K	NX665K	NX666K	NX667K	NX668K
Right	NX672K	NX673K	NX674K	NX675K	NX676K	NX677K	NX678K



### e.motion® UC Pro/PS Pro Tibia-Augments Lateral

4 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left + Right	NX612K	NX613K	NX614K	NX615K	NX616K	NX617K	NX618K

8 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX642K	NX643K	NX644K	NX645K	NX646K	NX647K	NX648K
Right	NX652K	NX653K	NX654K	NX655K	NX656K	NX657K	NX658K

12 mm							
Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX682K	NX683K	NX684K	NX685K	NX686K	NX687K	NX688K
Right	NX692K	NX693K	NX694K	NX695K	NX696K	NX697K	NX698K

# AESCULAP® e.motion® Pro SYSTEM

19 | IMPLANT MATRIX

## AS e.motion® Pro Implant Matrix – Femoral implants



### Femur FP/UC cemented

Types:	F2	F3	F4	F4N	F5	F5N	F6	F6N	F7	F8
Left	N0502Z	N0503Z	N0504Z	N0817Z	N0505Z	N0818Z	N0506Z	N0819Z	N0507Z	N0508Z
Right	N0602Z	N0603Z	N0604Z	N0917Z	N0605Z	N0918Z	N0606Z	N0919Z	N0607Z	N0608Z



### AS Femur PS Pro cemented

Types:	F2	F3	F4N	F4	F5N	F5	F6N	F6	F7	F8
Left	NX700Z	NX701Z	NX702Z	NX703Z	NX704Z	NX705Z	NX706Z	NX707Z	NX708Z	NX709Z
Right	NX750Z	NX751Z	NX752Z	NX753Z	NX754Z	NX755Z	NX756Z	NX757Z	NX758Z	NX759Z

## AS e.motion® Pro Implant Matrix – Tibial implants



### Tibia UC Pro/PS Pro cemented

Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX732Z	NX733Z	NX734Z	NX735Z	NX736Z	NX737Z	NX738Z
Right	NX782Z	NX783Z	NX784Z	NX785Z	NX786Z	NX787Z	NX788Z

## AS e.motion® Pro Implant Matrix – Tibial implants



### AS Rotation axis for meniscus components SW 4.5

#### Types:

Height/mm	10	12	14	16	18	20	22	24
UC	NR801Z	NR811Z	NR821Z	NR831Z	NR841Z	NR851Z	--	--
PS	NB800Z	NB810Z	NB820Z	NB830Z	NB840Z	NB850Z	NB860Z	NB870Z



### AS e.motion® UC Pro/PS Pro Tibia-Augments Medial

#### 4 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left + Right	NX602Z	NX603Z	NX604Z	NX605Z	NX606Z	NX607Z	NX608Z

#### 8 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX622Z	NX623Z	NX624Z	NX625Z	NX626Z	NX627Z	NX628Z
Right	NX632Z	NX633Z	NX634Z	NX635Z	NX636Z	NX637Z	NX638Z

#### 12 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX662Z	NX663Z	NX664Z	NX665Z	NX666Z	NX667Z	NX668Z
Right	NX672Z	NX673Z	NX674Z	NX675Z	NX676Z	NX677Z	NX678Z



### AS e.motion® UC Pro/PS Pro Tibia-Augments Lateral

#### 4 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left + Right	NX612Z	NX613Z	NX614Z	NX615Z	NX616Z	NX617Z	NX618Z

#### 8 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX642Z	NX643Z	NX644Z	NX645Z	NX646Z	NX647Z	NX648Z
Right	NX652Z	NX653Z	NX654Z	NX655Z	NX656Z	NX657Z	NX658Z

#### 12 mm

Types:	T2	T3	T4	T5	T6	T7	T8
Left	NX682Z	NX683Z	NX684Z	NX685Z	NX686Z	NX687Z	NX688Z
Right	NX692Z	NX693Z	NX694Z	NX695Z	NX696Z	NX697Z	NX698Z



# AESCULAP® e.motion® Pro SYSTEM

19 | IMPLANT MATRIX

## e.motion® Pro Implant Matrix – Meniscus Components



UC Pro – Left							
Types:	F2	F3	F4	F5	F6	F7	F8
10	NX402	NX403	NX404	NX405	NX406	NX407	NX408
12	NX412	NX413	NX414	NX415	NX416	NX417	NX418
14	NX422	NX423	NX424	NX425	NX426	NX427	NX428
16	NX432	NX433	NX434	NX435	NX436	NX437	NX438
18	NX442	NX443	NX444	NX445	NX446	NX447	NX448
20	NX452	NX453	NX454	NX455	NX456	NX457	NX458



PS Pro – Left							
	F2	F3	F4	F5	F6	F7	F8
	NX802	NX803	NX804	NX805	NX806	NX807	NX808
	NR812	NX813	NX814	NX815	NX816	NX817	NX818
	NX822	NX823	NX824	NX825	NX826	NX827	NX828
	NX832	NX833	NX834	NX835	NX836	NX837	NX838
	NX842	NX843	NX844	NX845	NX846	NX847	NX848
	NX852	NX853	NX854	NX855	NX856	NX857	NX858



UC Pro – Right							
Types:	F2	F3	F4	F5	F6	F7	F8
10	NX502	NX503	NX504	NX505	NX506	NX507	NX508
12	NX512	NX513	NX514	NX515	NX516	NX517	NX518
14	NX522	NX523	NX524	NX525	NX526	NX527	NX528
16	NX532	NX533	NX534	NX535	NX536	NX537	NX538
18	NX542	NX543	NX544	NX545	NX546	NX547	NX548
20	NX552	NX553	NX554	NX555	NX556	NX557	NX558



PS Pro – Right							
	F2	F3	F4	F5	F6	F7	F8
	NX902	NX903	NX904	NX905	NX906	NX907	NX908
	NX912	NX913	NX914	NX915	NX916	NX917	NX918
	NX922	NX923	NX924	NX925	NX926	NX927	NX928
	NX932	NX933	NX934	NX935	NX936	NX937	NX938
	NX942	NX943	NX944	NX945	NX946	NX947	NX948
	NX952	NX953	NX954	NX955	NX956	NX957	NX958

**NOTE:** All meniscus components include a rotation peg with AS coating.  
The AS rotation pegs can be combined with all AS and CoCr tibia components.

## e.motion® Pro Implant Matrix – Stems



### Tibia Extension Stem cemented

Types:	Ø 10 mm		Ø 12 mm		Ø 14 mm	
Length mm	52	92	52	92	52	92
Standard	NX060K	NX061K	NX062K	NX064K	NX063K	NX065K
AS	NX060Z	NX061Z	NX062Z	NX064Z	NX063Z	NX065Z



### Tibia Stems cementless

Types:	Ø 10 mm		Ø 12 mm		Ø 14 mm	
Length mm	92	132	92	132	92	132
Standard	NX082K	NX083K	NX084K	NX086K	NX085K	NX087K
AS	NX082Z	NX083Z	NX084Z	NX086Z	NX085Z	NX087Z

# AESCULAP® e.motion® Pro SYSTEM

20 | LITERATURE

- (1) Eiff W. Prozessoptimierung und Kostensenkung. HCM. 2016 Dec;7:34-7.
- (2) Grupp TM, Schroeder C, Kyun Kim T, Miehke RK, Fritz B5, Jansson V, Utzschneider S. Biotribology of a mobile bearing posterior stabilised knee design--effect of motion restraint on wear, tibio-femoral kinematics and particles. J Biomech. 2014 Jul 18;47(5):2415-23. Epub 2014 Apr 30.
- (3) Scheibel MT, Thomas M, von Salis-Soglio G. Operative Zugangswege in der Primärendoprothetik des Kniegelenks. Orthopäde. 2002;31:934-46. doi:10.1007/s00132-002-0383-0.
- (4) Amirfeyz R, Bannister G. The effect of bone porosity on the shear strength of the bone-cement interface. Int. Orthop. 2009 Jun;33(3):843-6.
- (5) Seeger JB1, Jaeger S, Bitsch RG, Mohr G, Rohner E, Clarius M. The effect of bone lavage on femoral cement penetration and interface temperature during Oxford unicompartmental knee arthroplasty with cement. J Bone Joint Surg Am. 2013 Jan 2;95(1):48-53.
- (6) Schlegel UJ1, Puschel K, Morlock MM, Nagel K. An in vitro comparison of tibial tray cementation using gun pressurization or pulsed lavage. 2014 May;38(5):967-71.
- (7) Norton MR, Eyres KS. Irrigation and suction technique to ensure reliable cement penetration for Total Knee Arthroplasty. J Arthroplasty. 2000 Jun;15(4):468-74.
- (8) British Orthopaedic Association and British Association for Surgery of the Knee. Knee Replacement: a guide to good practice: London: British Orthopaedic Association.
- (9) Vaninbrouckx M, Labey L, Innocenti B, Bellemans J. Cementing the femoral component in total knee arthroplasty: which technique is the best? Knee. 2009 Aug;16(4):265-8. doi: 10.1016/j.knee.2008.11.015.
- (10) De Baets T, Waelput W, Bellemans J. Analysis of third body particles generated during Total Knee Arthroplasty: is metal debris an issue? Knee. 2008 Mar;15(2):95-7. 2011.

NOTES

Lined area for taking notes, consisting of multiple horizontal dotted lines.

# AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | [www.aesculap.com](http://www.aesculap.com)

The main product trademark "AESCULAP" and the product trademarks "Aesculap OrthoTray", "Aesculap Reset", "e.motion", "Aesculap OrthoPilot" and "PLASMAPORE" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.