


Medibox® Sharps Disposal Containers

TECHNICAL DATA SHEET




Device variants	Medibox® 0.8 L	Medibox® 2.4 L	Medibox® 4.7 L	Medibox® 5.7 L	Medibox® 6.8 L	Medibox® 9.1 L
REF	9193424	9193600	9193618	9193627	9193639	9193640
Volume total	0.8 L	2.4 L	4.7 L	5.7 L	6.8 L	9.1 L
max. filling volume	0.67 L	2.0 L	3.9 L	4.8 L	5.7 L	7.7 L



Size	A: height (incl. lid)	18 cm	20 cm	20.5 cm	25.7 cm	50 cm	32.5 cm
	B: width (top side)	9.8 cm	15.1 cm	19 cm	19 cm	16.7 cm	21.5 cm
	B: depth (top side)	9.8 cm	15.1 cm	19 cm	19 cm	16.7 cm	21.5 cm
	C: base dimensions	ø 7.5 cm	11 x 11 cm	15 x 15 cm	13.9x13.9 cm	ø 11 cm	17 x 17 cm
	wall thickness	1.9 mm	2 mm	1.2 mm	1.2 mm	2 mm	1.9 mm
Pieces per carton		18 pieces	27 pieces	25 pieces	25 pieces	22 pieces	13 pieces
Carton dimensions L/W/H cm		39.5 x 29.5 x 16	79 x 39 x 18	79 x 39 x 21.8	79 x 39 x 21.8	79 x 39 x 33.5	79 x 39 x 33.5

Material	Polypropylene	
Colour can/cap	yellow/white	
Closure	cap for temporary closure final lock mechanism	
Certificates	UN type approval acc. ADR 2015 and 2008/68/EC	
	Kitemark approval acc. ISO 23907	
Product classification	consumer product sharps collectors are not classified as medical devices acc. to European Medical Device Regulation (EU) 2017/745	

- Disposal of
1. luer and luer lock needles
 2. pen needles
 3. various medical sharps



Medibox® 0.8 L – 9.1 L