

Actreen® Hi-Lite Set

... to use anywhere, anytime!

Patient benefits



DISCREET



LIGHTWEIGHT



CONVENIENT,
EASY TO USE AND
READY TO USE



SAFE



ENVIRONMENTALLY
FRIENDLY



Actreen® Hi-Lite Set

Ordering information

Tiemann 37 cm	CH	Diameter	REF	Unit/Box
	■ 10	ø 3,3 mm	242110*	30 + 1 discreet pouch
	□ 12	ø 4,0 mm	242112*	30 + 1 discreet pouch
	■ 14	ø 4,7 mm	242114*	30 + 1 discreet pouch
	■ 16	ø 5,3 mm	242116*	30 + 1 discreet pouch
	■ 18	ø 6,0 mm	242118*	30

Nelaton 37 cm	CH	Diameter	REF	Unit/Box
	■ 10	ø 3,3 mm	242210*	30 + 1 discreet pouch
	□ 12	ø 4,0 mm	242212*	30 + 1 discreet pouch
	■ 14	ø 4,7 mm	242214*	30 + 1 discreet pouch
	■ 16	ø 5,3 mm	242216*	30 + 1 discreet pouch
	■ 18	ø 6,0 mm	242218*	30

* For each country a specific letter is added at the end of the reference



The packaging of Actreen® Hi-Lite Set is composed of recyclable materials. The polyethylene plastic material turns into water vapour and carbon dioxide if incinerated. Discard the catheter according to the current local regulations.

Manufactured and distributed by:
B. Braun Medical | 204 avenue du Maréchal Juin | 92100 Boulogne-Billancourt | FRANCE
Tél. 01 41 10 53 00 | Fax 01 41 10 53 99 | www.bbraun.com
B. Braun Medical | SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun. Non binding documents and photographs.

Nr. ZG00744 - 01/2016



Actreen® Hi-Lite Set

ANYWHERE, ANYTIME
Self intermittent urinary set

Actreen® Hi-Lite Set

An invisible solution all day long...

A discreet and ready to use set...

...safe and convenient...



Everyone's daily life is filled with both personal and professional activities!

If you suffer from urinary disorders and are in need of a urinary set then you will be looking for a **discreet, simple and ready to use** product to use on a daily basis.

Actreen® Hi Lite Set is the latest generation of set (a catheter sealed to a bag) especially developed to make your life easier anywhere and anytime.

Whether you are in your car, at work or even at a restaurant, this set is ideal due to its **discreet and convenient** features.

We are all looking for ways to make our day easier regardless of the daily worries!



Safe material:
PVC Free, DEHP Free
Non traumatic eyes²

Simple to open:
Opening adapted to each patient's level of dexterity¹

Discreet, anywhere...¹
Actreen® Hi Lite Set can be easily folded⁴, and stored inside the provided pouch for more discretion throughout the day

LIGHTWEIGHT FOLDABLE

LOW DEXTERITY

GOOD DEXTERITY

Let yourself discover this unique solution which will allow you to live your life to the fullest.

¹ Report n° 7.5.AB.054.1 + Results report of the USABILITY TEST Product: Actreen® Hi Lite Set (nurses/patients)

² Report n° 7.5.AB.054.1

³ Results report of the USABILITY TEST Product: Actreen® Hi Lite Set (nurses/patients)

⁴ The Ch 18 cannot be folded, and no discreet pouch is included.

⁵ For safety and hygienic reasons, we advise you to use the product in the shortest possible time after opening.

⁶ Report n° 7.5.AC.059.2

⁷ Reports n° 7.5.AB.057.1, N° 7.5.AB.058.1, 7.5.AB.059.1, N° 7.5.AB.060.1