


# Folitetex<sup>®</sup> Sensitive tie-on

MOUTH AND NOSE COVERING SURGICAL MASK | TECHNICAL DATA



## PRODUCT CLASSIFICATION

Medical Device Class I

## QUALITY CERTIFICATES

ISO 9001, ISO 13485

## APPLIED STANDARDS

EN 14683, ISO 22609, ISO 10993

## INTENDED USE

ACC. EN 14683:2019

The main intended use of medical face masks is to protect the patient from infective acc. EN 14683:2019 agents and, additionally, in certain circumstances to protect the wearer against splashed of potentially contaminated liquids. Medical face masks may also be intended to be worn by patients and other persons to reduce the risk of spread of infections, particularly in epidemic or pandemic situations.

## MANUFACTURER



priMED Medical Products Inc.,  
200, 2003 91 St. SW, Edmonton, AB, Canada T6X 0W8







## IMPORTER



B. Braun Melsungen AG  
Carl-Braun-Str. 1, 34212 Melsungen  
Germany

# Folitetex<sup>®</sup> Sensitive tie-on

## MOUTH AND NOSE COVERING SURGICAL MASK | TECHNICAL DATA

<b>MEDICAL DEVICE INFORMATION</b>	Medical device class I, EN 14683	
<b>MASK DESIGN</b>	     	
<b>MATERIAL</b>	Colour Layers Fixation Outer layer Filter layer Inner layer Nose clip Ties	white 3ply, folded, ultrasonically welded ties, flexible nose clip spunbond polypropylene, white melt-blown polypropylene non-woven polyethylene/polypropylene bi-component, white aluminum strip spunbond polypropylene
<b>PHYSICAL PROPERTIES</b>	Free of natural rubber latex	
<b>DIMENSIONS</b>	Protection type Barrier efficiency Breathability Splash resistance Microbial cleanliness	type IIR ≥ 98 % bacterial filtration efficiency BFE < 60 Pa/cm <sup>2</sup> differential pressure for air permeation withstanding penetration of artificial blood at a pressure of ≥ 16 kPa ≤ 30 CFU/g (colony forming units)
<b>LOGISTIC INFORMATION</b>	Mask Ties (total) Nose Fixation Dispenser pack Transportation carton Storage conditions	180 mm width/93 mm height ≥ 91 cm length 140 mm length 6073312 50 pieces                      190 x 110 x 150 mm (L x W x H) 6 dispenser packs            460 x 230 x 20 mm (L x W x H) protect from direct heat, lighting or sunlight/UV radiation, store clean and dry at ambient temperature